

Package ‘angstromATE’

October 20, 2024

Type Package

Title Imports Log Files from Angstrom Engineering Thermal Evaporator

Version 0.1.3

Author Thomas Gredig [aut, cre, cph] (<<https://orcid.org/0000-0002-5824-7626>>)

Maintainer Thomas Gredig <tgredig@csulb.edu>

Description Opens and imports log files from Angstrom Engineering Thermal Evaporator and extracts basic characteristics, such as base pressure, time of the evaporation. It can visualize the deposition observables for review.

License GPL (>= 3)

BugReports <https://github.com/thomasgredig/angstromATE/issues>

Encoding UTF-8

RoxygenNote 7.3.2

Imports XML, stringr, utils

Suggests testthat (>= 3.0.0)

Config/testthat/edition 3

NeedsCompilation no

Repository CRAN

Date/Publication 2024-10-20 21:50:02 UTC

Contents

| | |
|---------------------------|---|
| ATE.complete | 2 |
| ATE.import | 2 |
| ATE.info | 3 |
| ATE.sampleFiles | 4 |
| ATE.status | 4 |
| conv2seconds | 5 |

| | |
|--------------|----------|
| Index | 6 |
|--------------|----------|

ATE.complete

ATE Deposition Summary

Description

Reads an XML Status file from the ATE thermal evaporator and returns the procedural timeline.

Usage

```
ATE.complete(filename, summaryOnly = FALSE)
```

Arguments

| | |
|-------------|--------------------------------------------------|
| filename | XML Status file from Angstrom Thermal Evaporator |
| summaryOnly | logical, if TRUE, returns summary for one layer |

Details

Extracts information about the deposition thickness from the completed status XML file at the end of the deposition.

Value

deposition thickness, rate, ramping times, and actions

Author(s)

Thomas Gredig

Examples

```
fileName <- ATE.sampleFiles('_Complete_')
ATE.complete(fileName, TRUE)
ATE.complete(fileName)
```

ATE.import*Imports Angstrom Engineering Thermal Evaporator Log Data*

Description

Imports Angstrom Engineering Thermal Evaporator Log Data

Usage

```
ATE.import(filename)
```

Arguments

filename CSV filename including path for the ATE log file

Value

data frame with around 50 variables and rows that represent time; the variables include Date, Time, SubstrateShutterOpen, ChamberPressure and many other parameters.

Author(s)

Thomas Gredig

Examples

```
fileName = ATE.sampleFiles("csv")[1]
d = ATE.import(fileName)
head(d)
```

ATE.info

Basic Information from Thermal Evaporator Deposition

Description

Imports data from the CSV log file of an Angstrom Engineering Thermal Evaporator. It extracts information during the deposition; i.e. while the shutter is open. It returns a condensed version of the deposition parameters.

Usage

```
ATE.info(filename, verbose = FALSE)
```

Arguments

filename full path of ATE Log file
verbose set to TRUE to get additional information

Value

list with information during the deposition that includes the thickness, the deposition time in seconds, the starting date, the substrate heater temperature, the material deposition temperature, tooling factor, base pressure, the pressure at start of the deposition, maximum pressure, and the material name

Author(s)

Thomas Gredig

See Also

[ATE.import()]

Examples

```
fileName = ATE.sampleFiles("csv")[1]
d = ATE.info(fileName,TRUE)
head(d)
```

| | |
|-----------------|-----------------------------|
| ATE.sampleFiles | <i>Sample ATE file list</i> |
|-----------------|-----------------------------|

Description

Returns a list of sample thermal evaporator log files, mostly for testing.

Usage

```
ATE.sampleFiles(filePattern = "*")
```

Arguments

filePattern pattern to limit the files

Value

list of sample data files with log information

Examples

```
ATE.sampleFiles()
ATE.sampleFiles('_Status')
```

| | |
|------------|--------------------------|
| ATE.status | <i>ATE Recipe Status</i> |
|------------|--------------------------|

Description

Reads an XML Status file from the ATE thermal evaporator and returns the procedural timeline.

Usage

```
ATE.status(filename)
```

Arguments

filename path and filename of XML Status file from Angstrom Thermal Evaporator

Value

data frame with description steps, start and end times

Author(s)

Thomas Gredig

Examples

```
fileName <- ATE.sampleFiles('_Status')
ATE.status(fileName)
```

conv2seconds *Convert Time String to Numeric*

Description

Convert Time String to Numeric

Usage

```
conv2seconds(strTime)
```

Arguments

strTime a string with time

Value

a numeric value in units of seconds

Author(s)

Thomas Gredig

Examples

```
conv2seconds("00:35:40.1816298")
conv2seconds("00:35:40.1816298") - conv2seconds("00:36:40.1816298")
conv2seconds("1.19:07:06.5180408")
```

Index

ATE.complete, [2](#)
ATE.import, [2](#)
ATE.info, [3](#)
ATE.sampleFiles, [4](#)
ATE.status, [4](#)

conv2seconds, [5](#)