Package 'ptairData'

August 8, 2024

Title PTR-TOF-MS volatolomics raw datasets from exhaled air and cell culture headspace

Version 1.13.0

Description The package ptairData contains two raw datasets from

Proton-Transfer-Reaction Time-of-Flight mass spectrometer acquisitions (PTR-TOF-MS), in the HDF5 format. One from the exhaled air of two volunteer healthy individuals with three replicates, and one from the cell culture headspace from two mycobacteria species and one control (culture medium only) with two replicates. Those datasets are used in the examples and in the vignette of the ptairMS package (PTR-TOF-MS data pre-processing). There are also used to gererate the ptrSet in the ptairMS data : exhaledPtrset and mycobacteriaSet

License GPL-3

Encoding UTF-8

RoxygenNote 7.1.1

BugReports https://github.com/camilleroquencourt/ptairData/issues

biocViews ExperimentData, MassSpectrometryData, CellCulture

Suggests knitr, rmarkdown, testthat, BiocStyle

VignetteBuilder knitr

Imports rhdf5, signal, stats, graphics

git_url https://git.bioconductor.org/packages/ptairData

git_branch devel

git_last_commit 97ff01c

git_last_commit_date 2024-04-30

Repository Bioconductor 3.20

Date/Publication 2024-08-08

Author camille Roquencourt [aut, cre]

Maintainer camille Roquencourt <camille.roquencourt@hotmail.fr>

Contents

Index

ptairData-package

ptairData: PTR-TOF-MS volatolomics raw datasets from exhaled air and cell culture headspace

Description

The package ptairData contains two raw datasets from Proton-Transfer-Reaction Time-of-Flight mass spectrometer acquisitions (PTR-TOF-MS), in the HDF5 format: one from the exhaled air of two volunteer healthy individuals, and one from the cell culture headspace from two mycobacteria species and one control (culture medium only). Those datasets are used in the examples and in the vignette from the ptairMS package (PTR-TOF-MS data pre-processing).

Author(s)

Maintainer: camille Roquencourt <camille.roquencourt@hotmail.fr>

3

Index

* internal
 ptairData-package, 2

ptairData(ptairData-package), 2
ptairData-package, 2